

# Recombinant Agrobacterium sp. aroA/EPSPS protein ,No Tag

|                     |  |
|---------------------|--|
| 货号                  | YP-Ab-17780  |
| 种属                  | Agrobacterium sp. (strain CP4)   |
| 纯度                  | >90% as determined by SDS-PAGE   |
| 体内研究                | Please contact with the lab for this information.  |
| 序列                  | A DNA sequence encoding the Agrobacterium sp. aroA (His3-Ala455) was fused without Tag   |
| Human Swiss Prot No | Q9R4E4   |
| 组成                  | Supplied as solution form in 50 mM Tris-HCl pH 7.5, 150 mM NaCl./ Lyophilized from 50 mM Tris-HCl pH 7.5, 150 mM NaCl.   |
| 来源                  | E.coli   |
| 稀释                  | Reconstitute in sterile water for a stock solution   |
| 浓度                  | >90% as determined by SDS-PAGE   |
| 储存                  | Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week . Store at -20 to -80 °C for twelve months from the date of receipt. |
| 其他名称                | EPSPS, EPSP synthase, 5-enolpyruvylshikimate-3-phosphate synthase  |
| 有关注意事项              | Avoid repeated freezing and thawing!   |
| 使用建议                | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.  |

## Products Images