



MIRO2 Rabbit pAb

Catalog No	YP-Ab-18449
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	RHOT2 ARHT2 C16orf39
Protein Name	Mitochondrial Rho GTPase 2 (MIRO-2) (hMiro-2) (Ras homolog gene family member T2)
Immunogen	Synthesized peptide derived from human MIRO2
Specificity	This antibody detects endogenous levels of MIRO2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	68kD
Cell Pathway	Mitochondrion outer membrane ; Single-pass type IV membrane protein . Colocalizes with MGARP and RHOT2 at the mitochondria.
Tissue Specificity	Mitochondrion outer membrane ; Single-pass type IV membrane protein . Colocalizes with MGARP and RHOT2 at the mitochondria.
Function	Mitochondrial GTPase involved in mitochondrial trafficking . Probably involved in control of anterograde transport of mitochondria and their subcellular distribution .
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

