



ALG6 Rabbit pAb

Catalog No	YP-Ab-18451
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	ALG6 My046
Protein Name	Dolichyl pyrophosphate Man9GlcNAc2 alpha-1,3-glucosyltransferase (Asparagine-linked glycosylation protein 6 homolog) (Dol-P-Glc:Man(9)GlcNAc(2)-PP-Dol alpha-1,3-glucosyltransferase) (Dolichyl-P-Glc:Man9GlcNAc2-PP-dolichyl glucosyltransferase)
Immunogen	Synthesized peptide derived from human ALG6
Specificity	This antibody detects endogenous levels of ALG6 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	56kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	Adds the first glucose residue to the lipid-linked oligosaccharide precursor for N-linked glycosylation. Transfers glucose from dolichyl phosphate glucose (Dol-P-Glc) onto the lipid-linked oligosaccharide Man(9)GlcNAc(2)-PP-Dol.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images