



Collagen XV α 1 (Cleaved-Asn1212) mouse mAb

Catalog No	YP-mAb-04380
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	COL15A1
Protein Name	Collagen XV α 1 (Cleaved-Asn1212)
Immunogen	Synthesized peptide derived from human Collagen XV α 1 (Cleaved-Asn1212)
Specificity	This antibody detects endogenous levels of Human,Mouse Collagen XV α 1 (Cleaved-Asn1212, protein was cleaved amino acid sequence between 1211-1212)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	Collagen alpha-1(XV) chain [Cleaved into: Endostatin (Endostatin-XV;Related to endostatin;Restin)]
Observed Band	20 150kD
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	Expressed predominantly in internal organs such as adrenal gland, pancreas and kidney.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images