



KLK1/2/3 (Cleaved-Ile25) mouse mAb

| | |
|----------------------------------|---|
| Catalog No | YP-mAb-04387 |
| Isotype | IgG |
| Reactivity | Human;Rat;Mouse; |
| Applications | WB |
| Gene Name | KLK1/2/3 |
| Protein Name | KLK1/2/3 (Cleaved-Ile25) |
| Immunogen | Synthesized peptide derived from human KLK1/2/3 (Cleaved-Ile25) |
| Specificity | This antibody detects endogenous levels of Human KLK1/2/3 (Cleaved-Ile25, protein was cleaved amino acid sequence between 24-25) |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Source | Monoclonal, Mouse,IgG |
| Purification | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution | WB 1:500-1:2000 |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| Synonyms | Kallikrein-2 (EC 3.4.21.35;Glandular kallikrein-1;hGK-1;Tissue kallikrein-2) |
| Observed Band | 25kD |
| Cell Pathway | |
| Tissue Specificity | |
| Function | |
| Background | |
| matters needing attention | Avoid repeated freezing and thawing! |
| Usage suggestions | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



Products Images