



# MCFD2 (Cleaved-Glu27) mouse mAb

<b>Catalog No</b>	YP-mAb-04395
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Mouse;Rat
<b>Applications</b>	WB
<b>Gene Name</b>	MCFD2 SDNSF
<b>Protein Name</b>	MCFD2 (Cleaved-Glu27)
<b>Immunogen</b>	Synthesized peptide derived from human MCFD2 (Cleaved-Glu27)
<b>Specificity</b>	This antibody detects endogenous levels of Human,Mouse,Rat MCFD2 (Cleaved-Glu27, protein was cleaved amino acid sequence between 26-27 )
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	Monoclonal, Mouse,IgG
<b>Purification</b>	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-1:2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	Multiple coagulation factor deficiency protein 2 (Neural stem cell-derived neuronal survival protein)
<b>Observed Band</b>	14 16kD
<b>Cell Pathway</b>	Endoplasmic reticulum-Golgi intermediate compartment . Endoplasmic reticulum . Golgi apparatus .
<b>Tissue Specificity</b>	
<b>Function</b>	
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images