



IIGP5 Monoclonal Antibody

Catalog No	YP-mAb-05102
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	IRGC IIGP5 IRGC1
Protein Name	Interferon-inducible GTPase 5 (EC 3.6.5.-) (Immunity-related GTPase cinema 1)
Immunogen	Synthesized peptide derived from human protein . at AA range: 350-430
Specificity	IIGP5 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	50kD
Cell Pathway	membrane,
Tissue Specificity	Brain,Testis,
Function	catalytic activity:GTP + H(2)O = GDP + phosphate.,similarity:Belongs to the interferon-inducible GTPase family.,
Background	catalytic activity:GTP + H(2)O = GDP + phosphate.,similarity:Belongs to the interferon-inducible GTPase family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images