



B3A2 Monoclonal Antibody

Catalog No	YP-mAb-05321
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SLC4A2 AE2 EPB3L1 HKB3 MPB3L
Protein Name	Anion exchange protein 2 (AE 2) (Anion exchanger 2) (Non-erythroid band 3-like protein) (BND3L) (Solute carrier family 4 member 2)
Immunogen	Synthesized peptide derived from human protein . at AA range: 830-910
Specificity	B3A2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	136kD
Cell Pathway	Membrane; Multi-pass membrane protein.
Tissue Specificity	Brain,
Function	alternative products:Additional isoforms seem to exist,function:Plasma membrane anion exchange protein of wide distribution.,similarity:Belongs to the anion exchanger (TC 2.A.31) family.,
Background	alternative products:Additional isoforms seem to exist,function:Plasma membrane anion exchange protein of wide distribution.,similarity:Belongs to the anion exchanger (TC 2.A.31) family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images