



ABCF3 Monoclonal Antibody

Catalog No	YP-mAb-05354
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ABCF3
Protein Name	ATP-binding cassette sub-family F member 3
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	ABCF3 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	77kD
Cell Pathway	cell-cell adherens junction,membrane,
Tissue Specificity	Cervix,Muscle,Placenta,Skin,Subthalamic nucleus,Vascular sm
Function	caution:Lacks transmembrane domains and is probably not involved in transport.,similarity:Belongs to the ABC transporter family. EF3 subfamily.,similarity:Contains 2 ABC transporter domains.,
Background	caution:Lacks transmembrane domains and is probably not involved in transport.,similarity:Belongs to the ABC transporter family. EF3 subfamily.,similarity:Contains 2 ABC transporter domains.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images