



# MYBPH Monoclonal Antibody

<b>Catalog No</b>	YP-mAb-05775
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human;Rat;Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	MYBPH
<b>Protein Name</b>	Myosin-binding protein H (MyBP-H) (H-protein)
<b>Immunogen</b>	Synthesized peptide derived from human protein . at AA range: 30-110
<b>Specificity</b>	MYBPH Monoclonal Antibody detects endogenous levels of protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Source</b>	Monoclonal, Mouse,IgG
<b>Purification</b>	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-1:2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	52kD
<b>Cell Pathway</b>	myosin filament,
<b>Tissue Specificity</b>	Mainly expressed in the skeletal muscle. Slightly expressed in the left atrium and arteria mammaria interna.
<b>Function</b>	function: Binds to myosin; probably involved in interaction with thick myofilaments in the A-band.,similarity: Belongs to the immunoglobulin superfamily. MyBP family.,similarity: Contains 2 fibronectin type-III domains.,similarity: Contains 2 Ig-like C2-type (immunoglobulin-like) domains.,tissue specificity: Skeletal muscle.,
<b>Background</b>	function: Binds to myosin; probably involved in interaction with thick myofilaments in the A-band.,similarity: Belongs to the immunoglobulin superfamily. MyBP family.,similarity: Contains 2 fibronectin type-III domains.,similarity: Contains 2 Ig-like C2-type (immunoglobulin-like) domains.,tissue specificity: Skeletal muscle.,
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images