



PKHG1 Monoclonal Antibody

Catalog No	YP-mAb-05938
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	PLEKHG1 KIAA1209
Protein Name	Pleckstrin homology domain-containing family G member 1
Immunogen	Synthesized peptide derived from human protein . at AA range: 780-860
Specificity	PKHG1 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	152kD
Cell Pathway	
Tissue Specificity	Brain,Epithelium,Lymph,Pooled,Testicle,
Function	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 DH (DBL-homology) domain.,similarity:Contains 1 PH domain.,
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 DH (DBL-homology) domain.,similarity:Contains 1 PH domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images