



DEN5A Monoclonal Antibody

Catalog No	YP-mAb-06005
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	DENND5A KIAA1091 RAB6IP1
Protein Name	DENN domain-containing protein 5A (Rab6-interacting protein 1) (Rab6IP1)
Immunogen	Synthesized peptide derived from human protein . at AA range: 650-730
Specificity	DEN5A Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	141kD
Cell Pathway	Golgi apparatus membrane .
Tissue Specificity	Brain,Placenta,Skin,Uterus,
Function	function:May be involved in Rab6-mediated GTPase signaling.,similarity:Belongs to the RAB6IP1 family.,similarity:Contains 1 dDENN domain.,similarity:Contains 1 DENN domain.,similarity:Contains 1 PLAT domain.,similarity:Contains 1 uDENN domain.,similarity:Contains 2 RUN domains.,subunit:Probably interacts with RAB6A bound to GTP.,
Background	function:May be involved in Rab6-mediated GTPase signaling.,similarity:Belongs to the RAB6IP1 family.,similarity:Contains 1 dDENN domain.,similarity:Contains 1 DENN domain.,similarity:Contains 1 PLAT domain.,similarity:Contains 1 uDENN domain.,similarity:Contains 2 RUN domains.,subunit:Probably interacts with RAB6A bound to GTP.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images