



RIOK2 Monoclonal Antibody

Catalog No	YP-mAb-06137
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	RIOK2
Protein Name	Serine/threonine-protein kinase RIO2 (EC 2.7.11.1) (RIO kinase 2)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	RIOK2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	60kD
Cell Pathway	Cytoplasm . Exported out of the nucleus via its NES in a XPO1-dependent manner. .
Tissue Specificity	Epithelium,Placenta,
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. RIO-type Ser/Thr kinase family.,similarity:Contains 1 protein kinase domain.,
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. RIO-type Ser/Thr kinase family.,similarity:Contains 1 protein kinase domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images