



ZBT10 Monoclonal Antibody

Catalog No	YP-mAb-06403
Isotype	IgG
Reactivity	Human;Rat
Applications	WB
Gene Name	ZBTB10 RINZF RINZFC
Protein Name	Zinc finger and BTB domain-containing protein 10 (Zinc finger protein RIN ZF)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ZBT10 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	95kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,PCR rescued clones,Prostate,Testis,
Function	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 2 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 2 C2H2-type zinc fingers.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images