



# ZN786 Monoclonal Antibody

|                                  |                                                                                                                                                                                                                     |
|----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                | YP-mAb-06535                                                                                                                                                                                                        |
| <b>Isotype</b>                   | IgG                                                                                                                                                                                                                 |
| <b>Reactivity</b>                | Human;Rat;Mouse;                                                                                                                                                                                                    |
| <b>Applications</b>              | WB                                                                                                                                                                                                                  |
| <b>Gene Name</b>                 | ZNF786                                                                                                                                                                                                              |
| <b>Protein Name</b>              | Zinc finger protein 786                                                                                                                                                                                             |
| <b>Immunogen</b>                 | Synthesized peptide derived from part region of human protein                                                                                                                                                       |
| <b>Specificity</b>               | ZN786 Monoclonal Antibody detects endogenous levels of protein.                                                                                                                                                     |
| <b>Formulation</b>               | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.                                                                                                                                                      |
| <b>Source</b>                    | Monoclonal, Mouse,IgG                                                                                                                                                                                               |
| <b>Purification</b>              | The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                                |
| <b>Dilution</b>                  | WB 1:500-1:2000                                                                                                                                                                                                     |
| <b>Concentration</b>             | 1 mg/ml                                                                                                                                                                                                             |
| <b>Purity</b>                    | ≥90%                                                                                                                                                                                                                |
| <b>Storage Stability</b>         | -20°C/1 year                                                                                                                                                                                                        |
| <b>Synonyms</b>                  |                                                                                                                                                                                                                     |
| <b>Observed Band</b>             | 86kD                                                                                                                                                                                                                |
| <b>Cell Pathway</b>              | Nucleus .                                                                                                                                                                                                           |
| <b>Tissue Specificity</b>        | Cervix,Melanoma,PCR rescued clones,                                                                                                                                                                                 |
| <b>Function</b>                  | function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 16 C2H2-type zinc fingers., |
| <b>Background</b>                | function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 16 C2H2-type zinc fingers., |
| <b>matters needing attention</b> | Avoid repeated freezing and thawing!                                                                                                                                                                                |
| <b>Usage suggestions</b>         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.                                                                                   |



## Products Images