



ZNF71 Monoclonal Antibody

Catalog No	YP-mAb-06578
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF71 EZFIT
Protein Name	Endothelial zinc finger protein induced by tumor necrosis factor alpha (Zinc finger protein 71)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ZNF71 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	53kD
Cell Pathway	Nucleus .
Tissue Specificity	Highly expressed in placenta, followed by brain, testis, pancreas, heart, small intestine, muscle, uterus, prostate and peripheral blood leukocytes. Not detected in liver, lung, colon, stomach, salivary and thyroid gland.
Function	function:May be involved in transcriptional regulation.,induction:By tumor necrosis factor alpha (TNFalpha).,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 13 C2H2-type zinc fingers.,tissue specificity:Highly expressed in placenta, followed by brain, testis, pancreas, heart, small intestine, muscle, uterus, prostate and peripheral blood leukocytes. Not detected in liver, lung, colon, stomach, salivary and thyroid gland.,
Background	function:May be involved in transcriptional regulation.,induction:By tumor necrosis factor alpha (TNFalpha).,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 13 C2H2-type zinc fingers.,tissue specificity:Highly expressed in placenta, followed by brain, testis, pancreas, heart, small intestine, muscle, uterus, prostate and peripheral blood leukocytes. Not detected in liver, lung, colon, stomach, salivary and thyroid gland.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images