



RGMB Monoclonal Antibody

Catalog No	YP-mAb-06800
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	RGMB
Protein Name	RGM domain family member B (DRG11-responsive axonal guidance and outgrowth of neurite) (DRAGON)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	RGMB Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	48kD
Cell Pathway	Cell membrane ; Lipid-anchor, GPI-anchor . Membrane raft .
Tissue Specificity	Placenta,Testis,
Function	similarity:Belongs to the repulsive guidance molecule (RGM) family.,
Background	RGMB is a glycosylphosphatidylinositol (GPI)-anchored member of the repulsive guidance molecule family (see RGMA, MIM 607362) and contributes to the patterning of the developing nervous system (Samad et al., 2005 [PubMed 15671031]).[supplied by OMIM, Apr 2009],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images