



H6ST2 Monoclonal Antibody

Catalog No	YP-mAb-07242
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	HS6ST2 PSEC0092
Protein Name	Heparan-sulfate 6-O-sulfotransferase 2 (HS6ST-2) (EC 2.8.2.-)
Immunogen	Synthesized peptide derived from human protein . at AA range: 140-220
Specificity	H6ST2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	66kD
Cell Pathway	Membrane ; Single-pass type II membrane protein .
Tissue Specificity	Amygdala,Brain,Placenta,
Function	catalytic activity:3'-phosphoadenylyl sulfate + [heparan sulfate]-glucosamine = adenosine 3',5'-bisphosphate + [heparan sulfate]-glucosamine 6-sulfate.,function:6-O-sulfation enzyme which catalyzes the transfer of sulfate from 3'-phosphoadenosine 5'-phosphosulfate (PAPS) to position 6 of the N-sulfoglucosamine residue (GlcNS) of heparan sulfate.,similarity:Belongs to the sulfotransferase 6 family.,
Background	Heparan sulfate proteoglycans are ubiquitous components of the cell surface, extracellular matrix, and basement membranes, and interact with various ligands to influence cell growth, differentiation, adhesion, and migration. This gene encodes a member of the heparan sulfate (HS) sulfotransferase gene family, which catalyze the transfer of sulfate to HS. Different family members and isoforms are thought to synthesize heparan sulfates with tissue-specific structures and functions. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images