



CC14C Monoclonal Antibody

Catalog No	YP-mAb-07250
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CDC14C CDC14B2
Protein Name	Dual specificity protein phosphatase CDC14C (EC 3.1.3.16) (EC 3.1.3.48) (CDC14 cell division cycle 14 homolog C)
Immunogen	Synthesized peptide derived from human protein . at AA range: 450-530
Specificity	CC14C Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	60kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass membrane protein . Retains its endoplasmic reticulum localization during mitosis. .
Tissue Specificity	
Function	
Background	Dual-specificity phosphatase. Preferentially dephosphorylates proteins modified by proline-directed kinases (By similarity).
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images