



AF1L2 Monoclonal Antibody

Catalog No	YP-mAb-07688
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	AFAP1L2 KIAA1914 XB130
Protein Name	Actin filament-associated protein 1-like 2 (AFAP1-like protein 2)
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	AF1L2 Monoclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	89kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Detected in spleen and thyroid, and at lower levels in kidney, brain, lung and pancreas.
Function	function:May play a role in a signaling cascade by enhancing the kinase activity of SRC. Contributes to SRC-regulated transcription activation.,PTM:Tyrosine phosphorylated (by SRC).,similarity:Contains 2 PH domains.,subunit:Interacts with SRC. Interacts with LCK when tyrosine phosphorylated.,tissue specificity:Detected in spleen and thyroid, and at lower levels in kidney, brain, lung and pancreas.,
Background	function:May play a role in a signaling cascade by enhancing the kinase activity of SRC. Contributes to SRC-regulated transcription activation.,PTM:Tyrosine phosphorylated (by SRC).,similarity:Contains 2 PH domains.,subunit:Interacts with SRC. Interacts with LCK when tyrosine phosphorylated.,tissue specificity:Detected in spleen and thyroid, and at lower levels in kidney, brain, lung and pancreas.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images