



DJC24 mouse mAb

Catalog No	YP-mAb-09188
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	DNAJC24 DPH4 ZCSL3
Protein Name	DnaJ homolog subfamily C member 24 (CSL-type zinc finger-containing protein 3) (DPH4 homolog)
Immunogen	Synthesized peptide derived from human DJC24 AA range: 39-89
Specificity	This antibody detects endogenous levels of human DJC24
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm, cytoskeleton .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images