



EPMIP mouse mAb

Catalog No	YP-mAb-09192
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	EPM2AIP1 KIAA0766 My007
Protein Name	EPM2A-interacting protein 1 (Laforin-interacting protein)
Immunogen	Synthesized peptide derived from human EPMIP AA range: 180-230
Specificity	This antibody detects endogenous levels of human EPMIP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endoplasmic reticulum .
Tissue Specificity	Expressed in heart, brain, placenta, liver, pancreas, kidney and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images