



WBP4 mouse mAb

Catalog No	YP-mAb-09193
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	WBP4 FBP21 FBNP21
Protein Name	WW domain-binding protein 4 (WBP-4) (Formin-binding protein 21) (WW domain-containing-binding protein 4)
Immunogen	Synthesized peptide derived from human WBP4 AA range: 235-285
Specificity	This antibody detects endogenous levels of human WBP4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Nucleus speckle .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images