



TNPO3 mouse mAb

Catalog No	YP-mAb-09200
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TNPO3 IPO12
Protein Name	Transportin-3 (Importin-12) (Imp12) (Transportin-SR) (TRN-SR)
Immunogen	Synthesized peptide derived from human TNPO3 AA range: 601-651
Specificity	This antibody detects endogenous levels of human TNPO3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus envelope . Cytoplasm . Localizes to the nuclear envelope and annulate lamellae, which consists in stacks of endoplasmic reticulum membranes containing a high density of nuclear pores. .
Tissue Specificity	Expressed in skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images