



KRA15 mouse mAb

Catalog No	YP-mAb-09219
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	KRTAP1-5 KAP1.5 KRTAP1.5
Protein Name	Keratin-associated protein 1-5 (High sulfur keratin-associated protein 1.5) (Keratin-associated protein 1.5)
Immunogen	Synthesized peptide derived from human KRA15 AA range: 91-141
Specificity	This antibody detects endogenous levels of human KRA15
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	Expressed in the middle/upper portions of the hair cortex, in the region termed the keratogenous zone.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images