



TRUB1 mouse mAb

Catalog No	YP-mAb-09220
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRUB1 PUS4
Protein Name	Probable tRNA pseudouridine synthase 1 (EC 5.4.99.-)
Immunogen	Synthesized peptide derived from human TRUB1 AA range: 116-166
Specificity	This antibody detects endogenous levels of human TRUB1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Cytoplasm, cytosol . Catalyzes pseudouridylation of mRNAs in the nucleus. .
Tissue Specificity	Highly expressed in heart, skeletal muscle and liver. Expressed at lower levels in lung, small intestine, kidney and spleen.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images