



SEM3D mouse mAb

Catalog No	YP-mAb-09236
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SEMA3D UNQ760/PRO1491
Protein Name	Semaphorin-3D
Immunogen	Synthesized peptide derived from human SEM3D AA range: 442-492
Specificity	This antibody detects endogenous levels of human SEM3D
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images