



NOE1 mouse mAb

Catalog No	YP-mAb-09239
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	OLFM1 NOE1 NOEL1
Protein Name	Noelin (Neuronal olfactomedin-related ER localized protein) (Olfactomedin-1)
Immunogen	Synthesized peptide derived from human NOE1 AA range: 414-464
Specificity	This antibody detects endogenous levels of human NOE1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted . Cell junction, synapse . Endoplasmic reticulum . Cell projection, axon . Perikaryon .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images