



YIF1A mouse mAb

Catalog No	YP-mAb-09254
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	YIF1A 54TM HYIF1P YIF1
Protein Name	Protein YIF1A (54TMp) (YIP1-interacting factor homolog A)
Immunogen	Synthesized peptide derived from human YIF1A AA range: 101-151
Specificity	This antibody detects endogenous levels of human YIF1A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Golgi apparatus membrane ; Multi-pass membrane protein . Endoplasmic reticulum-Golgi intermediate compartment membrane ; Multi-pass membrane protein . Cycles between the endoplasmic reticulum and the endoplasmic reticulum-Golgi intermediate compartment. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images