



REBL1 mouse mAb

Catalog No	YP-mAb-09262
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RHEBL1
Protein Name	GTPase RhebL1 (Ras homolog enriched in brain like-1 c) (RhebL1c) (Ras homolog enriched in brain-like protein 1) (Rheb-like protein 1) (Rheb2)
Immunogen	Synthesized peptide derived from human REBL1 AA range: 32-82
Specificity	This antibody detects endogenous levels of human REBL1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Endomembrane system ; Lipid-anchor ; Cytoplasmic side . Cytoplasm .
Tissue Specificity	Ubiquitously expressed. Expression increased at least 2-fold in several tumor cell lines.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images