



CAND2 mouse mAb

Catalog No	YP-mAb-09268
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CAND2 KIAA0667 TIP120B
Protein Name	Cullin-associated NEDD8-dissociated protein 2 (Cullin-associated and neddylation-dissociated protein 2) (Epididymis tissue protein Li 169) (TBP-interacting protein of 120 kDa B) (TBP-interacting prote
Immunogen	Synthesized peptide derived from human CAND2 AA range: 291-341
Specificity	This antibody detects endogenous levels of human CAND2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in epididymis (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images