



RHG40 mouse mAb

Catalog No	YP-mAb-09275
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ARHGAP40 C20orf95
Protein Name	Rho GTPase-activating protein 40 (Rho-type GTPase-activating protein 40)
Immunogen	Synthesized peptide derived from human RHG40 AA range: 402-452
Specificity	This antibody detects endogenous levels of human RHG40
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images