



RD21L mouse mAb

Catalog No	YP-mAb-09297
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RAD21L1 RAD21L
Protein Name	Double-strand-break repair protein rad21-like protein 1
Immunogen	Synthesized peptide derived from human RD21L AA range: 501-551
Specificity	This antibody detects endogenous levels of human RD21L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Chromosome . In meiotic chromosomes, localized along axial elements in early meiosis: detectable on the axial elements in leptotene, and stays on the axial/lateral elements until mid pachytene. It then disappears and is replaced with RAD21. Compared to REC8, has mutually exclusive loading sites on the chromosomes: REC8 and RAD21L form distinct cohesin-enriched domains along the axial elements. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images