



H2A2C mouse mAb

Catalog No	YP-mAb-09301
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	HIST2H2AC H2AFQ
Protein Name	Histone H2A type 2-C (Histone H2A-GL101) (Histone H2A/q)
Immunogen	Synthesized peptide derived from human H2A2C AA range: 70-120
Specificity	This antibody detects endogenous levels of human H2A2C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus. Chromosome.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images