



HDHD1 mouse mAb

Catalog No	YP-mAb-09308
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	HDHD1 DXF68S1E FAM16AX GS1 HDHD1A
Protein Name	Pseudouridine-5'-monophosphatase (5'-PsiMPase) (EC 3.1.3.n6) (Haloacid dehalogenase-like hydrolase domain-containing protein 1) (Haloacid dehalogenase-like hydrolase domain-containing protein 1A) (Pro
Immunogen	Synthesized peptide derived from human HDHD1 AA range: 89-139
Specificity	This antibody detects endogenous levels of human HDHD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images