



PSD3 mouse mAb

Catalog No	YP-mAb-09310
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PSD3 EFA6R HCA67 KIAA0942
Protein Name	PH and SEC7 domain-containing protein 3 (Epididymis tissue protein Li 20mP) (Exchange factor for ADP-ribosylation factor guanine nucleotide factor 6) (Hepatocellular carcinoma-associated antigen 67) (
Immunogen	Synthesized peptide derived from human PSD3 AA range: 116-166
Specificity	This antibody detects endogenous levels of human PSD3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	115kD
Cell Pathway	Cell membrane . Cell projection, ruffle membrane . Cell junction, synapse, postsynaptic density . In interphase associated with the plasma membrane, in particular with membrane ruffling regions. .
Tissue Specificity	Isoform 2 is expressed in epididymis (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images