



SETD5 mouse mAb

Catalog No	YP-mAb-09321
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SETD5 KIAA1757
Protein Name	SET domain-containing protein 5
Immunogen	Synthesized peptide derived from human SETD5 AA range: 634-684
Specificity	This antibody detects endogenous levels of human SETD5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Chromosome . Localizes to active transcribed genes. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images