



CCD34 mouse mAb

Catalog No	YP-mAb-09325
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CCDC34 RAMA3 L15
Protein Name	Coiled-coil domain-containing protein 34 (Renal carcinoma antigen NY-REN-41)
Immunogen	Synthesized peptide derived from human CCD34 AA range: 249-299
Specificity	This antibody detects endogenous levels of human CCD34
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images