



MYCT1 mouse mAb

Catalog No	YP-mAb-09328
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MYCT1 MTLC MTMC1
Protein Name	Myc target protein 1 (Myc target in myeloid cells protein 1)
Immunogen	Synthesized peptide derived from human MYCT1 AA range: 114-164
Specificity	This antibody detects endogenous levels of human MYCT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Nucleus . Expressed in nuclei of hepatocellular carcinoma cell line BEL-7402 cells.
Tissue Specificity	Down-regulated in gastric cancer tissues.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images