



CJ107 mouse mAb

Catalog No	YP-mAb-09331
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	C10orf107
Protein Name	Uncharacterized protein C10orf107
Immunogen	Synthesized peptide derived from human CJ107 AA range: 7-57
Specificity	This antibody detects endogenous levels of human CJ107
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Cell projection, cilium, flagellum . Colocalized with pericentrin at centrosome of spermatocytes and round spermatids. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images