



NECP1 mouse mAb

Catalog No	YP-mAb-09335
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	NECAP1
Protein Name	Adaptin ear-binding coat-associated protein 1 (NECAP endocytosis-associated protein 1) (NECAP-1)
Immunogen	Synthesized peptide derived from human NECAP1 AA range: 109-159
Specificity	This antibody detects endogenous levels of human NECAP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasmic vesicle, clathrin-coated vesicle membrane . Cell membrane . Colocalizes with AP-2 at the plasma membrane. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images