



SPAT1 mouse mAb

Catalog No	YP-mAb-09349
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SPATA1
Protein Name	Spermatogenesis-associated protein 1 (Sperm-specific protein SP-2)
Immunogen	Synthesized peptide derived from human SPAT1 AA range: 369-419
Specificity	This antibody detects endogenous levels of human SPAT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasmic vesicle, secretory vesicle, acrosome .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images