



HERC2 mouse mAb

Catalog No	YP-mAb-09357
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	HERC2
Protein Name	E3 ubiquitin-protein ligase HERC2 (EC 6.3.2.-) (HECT domain and RCC1-like domain-containing protein 2)
Immunogen	Synthesized peptide derived from human HERC2 AA range: 646-696
Specificity	This antibody detects endogenous levels of human HERC2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	532kD
Cell Pathway	Cytoplasm. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole. Nucleus. Recruited to sites of DNA damage in response to ionizing radiation (IR) via its interaction with RNF8. May loose association with centrosomes during mitosis.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images