



CBLC mouse mAb

Catalog No	YP-mAb-09358
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CBLC CBL3 RNF57
Protein Name	Signal transduction protein CBL-C (RING finger protein 57) (SH3-binding protein CBL-3) (SH3-binding protein CBL-C)
Immunogen	Synthesized peptide derived from human CBLC AA range: 219-269
Specificity	This antibody detects endogenous levels of human CBLC
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	Ubiquitous.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images