



EX3L2 mouse mAb

Catalog No	YP-mAb-09360
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	EXOC3L2 XTP7
Protein Name	Exocyst complex component 3-like protein 2 (HBV X-transactivated gene 7 protein) (HBV XAg-transactivated protein 7)
Immunogen	Synthesized peptide derived from human EX3L2 AA range: 271-321
Specificity	This antibody detects endogenous levels of human EX3L2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images