



SIRPD mouse mAb

Catalog No	YP-mAb-09361
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SIRPD PTPNS1L2
Protein Name	Signal-regulatory protein delta (SIRP-delta) (Protein tyrosine phosphatase non-receptor type substrate 1-like 2)
Immunogen	Synthesized peptide derived from human SIRPD AA range: 105-155
Specificity	This antibody detects endogenous levels of human SIRPD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images