



HMHB1 mouse mAb

Catalog No	YP-mAb-09376
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	HMHB1
Protein Name	Minor histocompatibility protein HB-1 [Cleaved into: Minor histocompatibility antigen HB-1 (mHag HB-1)]
Immunogen	Synthesized peptide derived from human HMHB1 AA range: 1-24
Specificity	This antibody detects endogenous levels of human HMHB1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	Expressed in acute lymphoblastic leukemia B-cells and Epstein-Barr virus-transformed B-cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images