



PP6R1 mouse mAb

Catalog No	YP-mAb-09380
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PPP6R1 KIAA1115 PP6R1 SAPS1
Protein Name	Serine/threonine-protein phosphatase 6 regulatory subunit 1 (SAPS domain family member 1)
Immunogen	Synthesized peptide derived from human PP6R1 AA range: 136-186
Specificity	This antibody detects endogenous levels of human PP6R1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cytoplasm .
Tissue Specificity	Ubiquitous with higher expression in testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images