



GBRP mouse mAb

Catalog No	YP-mAb-09384
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	GABRP
Protein Name	Gamma-aminobutyric acid receptor subunit pi (GABA(A) receptor subunit pi)
Immunogen	Synthesized peptide derived from human GBRP AA range: 181-231
Specificity	This antibody detects endogenous levels of human GBRP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Cell junction, synapse, postsynaptic cell membrane; Multi-pass membrane protein. Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Most abundant in the uterus, also expressed in lung, thymus and prostate.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images