



TNR6C mouse mAb

Catalog No	YP-mAb-09391
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TNRC6C KIAA1582
Protein Name	Trinucleotide repeat-containing gene 6C protein
Immunogen	Synthesized peptide derived from human TNR6C AA range: 784-834
Specificity	This antibody detects endogenous levels of human TNR6C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	186kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images